



## Primer removal after a PCR reaction with Vivacon 2



PCR (Polymerase Chain Reaction) is one of the most versatile methods used today in molecular biology for a multitude of applications like preparing fragments for cloning or amplifying DNA sequences, for e.g. forensic analysis.

These applications may be more or less sensitive to the remaining components of the PCR reaction mixture. PCR reaction mixtures contain a variety of salts, free nucleotides, glycerol, proteins, and primers. Certain restriction enzymes as well as DNA ligase are particularly sensitive to the presence of contaminants in DNA samples. Due to this, most downstream applications will require some sort of PCR cleanup.

The PCR cleanup can be performed in a variety of ways:

1. Precipitation, using chemical solubility properties to selectively separate DNA from other sample components is a frequently used method. Its main disadvantage is the rather lengthy procedure along with the incomplete removal of co-precipitating buffer components and contaminants.
2. Chromatography, using size exclusion particles or affinity to glass to purify DNA from the PCR mixture components is a very sensitive and effective method. This technique is however costly, generally requires significant handling, and dilutes samples, so that they must be concentrated after elution from the matrix.
3. Ultrafiltration involves the isolation and concentration of PCR products using membrane devices. It is rapid, requires very little handling, achieves high recoveries, leaves DNA undamaged, and the concentrated DNA is free of contaminants that may inhibit downstream reactions.

Here, we show in an experiment the effective removal of primers using Vivacon 2 ultrafiltration devices and show that the 30 kDa Hydrosart membrane Vivacon 2 is effective at retaining 300 bp DNA fragments, while removing the 24 bp primers.

### Materials and Methods

In order to evaluate the effectiveness of primer removal after a PCR reaction, two PCR sample mock ups with excess of 24 bp primers (for optimal visibility on a SDS gel) were prepared in the following way: 20  $\mu\text{l}$  of a 75  $\mu\text{g}/\mu\text{l}$  24 bp primer solution in deionized water was added to 90  $\mu\text{l}$  of a 40  $\mu\text{g}/\mu\text{l}$  300 bp DNA fragment in deionized water and diluted with 320  $\mu\text{l}$  TE buffer, pH 8 to a final volume of 430  $\mu\text{l}$ .


200  $\mu\text{l}$  of a 75  $\mu\text{g}/\mu\text{l}$  24 bp primer solution in deionized water was added to 90  $\mu\text{l}$  of a 51  $\mu\text{g}/\mu\text{l}$  400 bp DNA fragment in deionized water and diluted with 320  $\mu\text{l}$  TE buffer, pH 8 to a final volume of 430  $\mu\text{l}$ .

For primer removal, 200  $\mu\text{l}$  of each mixture was diluted with 1800  $\mu\text{l}$  TE buffer, pH 8 in a Vivacon 2 with a 30 kDa membrane, then centrifuged for 20 minutes at 2500  $\times$  g. The sample was filled up again to 2 ml with TE buffer, pH 8 and centrifuged for additional 20 minutes at 2500  $\times$  g. The effectiveness of primer removal was analysed using a 12% TBE-Polyacrylamid SDS gel.

5  $\mu\text{l}$  samples of the initial sample, concentrate and concentrate after wash step were applied to the SDS gel. Duplicates were prepared of each step.

### Results

The SDS gel shows the effectiveness of primer removal with a 30 kDa Vivacon 2, with quantitative recoveries of the 300 bp and 400 bp PCR fragment in a 40 minute procedure. Using a 30 kDa Vivacon 2 device, primers and PCR reaction components can effectively be removed from a PCR sample containing 300 bp DNA fragments and larger for subsequent applications.



12 % TBE Polyacrylamid SDS gel  
30 kDa Vivacon 2

- Lane 1 300 bp DNA fragment + 24 bp Primer – original sample
- Lane 2 300 bp DNA fragment + 24 bp Primer – concentrate (1)
- Lane 3 300 bp DNA fragment + 24 bp Primer – concentrate (2)
- Lane 4 300 bp DNA fragment + 24 bp Primer – concentrate after wash (1)
- Lane 5 300 bp DNA fragment + 24 bp Primer – concentrate after wash (2)
- Lane 6 400 bp DNA fragment + 24 bp Primer – concentrate (1)
- Lane 7 400 bp DNA fragment + 24 bp Primer – concentrate (2)
- Lane 8 400 bp DNA fragment + 24 bp Primer – concentrate after wash (1)
- Lane 9 400 bp DNA fragment + 24 bp Primer – concentrate after wash (2)
- Lane 10 400 bp DNA fragment + 24 bp Primer – original sample

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

Phone +49.551.308.0  
Fax +49.551.308.3289  
www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178  
UK +44.1372.737159  
France +33.442.845600  
Italy +39.055.63.40.41  
Spain +34.91.3586102  
Japan +81.3.3740.5407

Specifications subject to change  
without notice. Printed and copyrighted  
by Sartorius Stedim Biotech GmbH  
W · G

Publication No.: SL-4047-e08111  
Order No.: 85032-535-16