Centrisart® I

Gebrauchsanleitung
 13229-E (5.000 D)
 13239-E (10.000 D)
 13249-E (20.000 D)
 13269-E (100.000 D)
 13279-E (300.000 D)

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
1. Remove the closing cap, allow the floater (inner tube) to slide out and stand it on a clean surface with open end down, avoiding contact with the membrane. Pour the sample (0.1–2.5 ml) into the centrifuge tube (outer tube, Fig. 1).

2. Re-insert the floater, membrane side down, and let it slide down onto the surface of the sample. Allow the centrifuge tube to stand (Fig. 2) for approx. 5 minutes to ensure complete wetting of the membrane. Exception: to filter noncoagulated blood, precentrifuge the tube for about 2 min. before inserting floater in order to avoid hemolysis. Do **not** re-insert closing cap for centrifugation.

3. Centrifuge in a laboratory centrifuge accepting up to 17 × 100 mm tubes (Fig. 3). Maximum relative centrifugal force: 2,500 g with swing head rotors, 2,000 g with angular head rotors (see graph). When filtering concentrated protein solutions, for example serum or blood, start with 5 minutes centrifugation at half of the maximum rpm. When filtering solutions containing lipids such as milk, only use an angular head rotor.

Maximum rpm for each rotor radius (r) between 10 and 20 cm.
Upper curve A:
for swing head rotors (2,500 g)
Lower curve B:
for angular head rotors (2,000 g)

4. Pipet out the filtrate (Fig. 4) as soon as possible after completion of centrifugation (filtrate will slowly diffuse back into the concentrate) ...

5. ... or remove the floater with opened forceps to have access to the concentrate (Fig. 5).

Note:
Not for use with organic solvents.

Sartorius ultrafilters contain glycerol. Due to its hygroscopicity, transparent spots of darker appearance may develop on the filters at higher levels of humidity. Such spots, however, do not have any adverse effects on the retention rate or on the flow rate of these filters.

If removal of glycerol (approx. 2 mg) is required, it can be effected by a 5-minute dip in water or other suitable liquid, or by precentrifugation, immediately prior to use. Do not allow membrane to dry out. Replace closing cap if losses by evaporation are to be avoided.

For the purpose of product development we reserve the right to make modifications.

Sartorius AG
Weender Landstrasse 94 - 108
37075 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

Specifications subject to change without notice.
Printed and copyrighted by Sartorius W101.46
Publication No.: SE-6002-a04057